



# DEPARAFFINIZATION STAND

Uniquely designed manual staining station is ideal for high throughput staining(s) like H&E, PAP Stain, Dehydration & Rehydration and all Histology & Cytology Stains. The jars of the staining station are resistant to xylene and alcohols.

Due to the narrowed mouth of jars, it saves up to 60% of reagent usage when compared to conventional glass jars. These high-quality solvent resistant jars hold approx. 250-275 ml of reagent. The station comes with 4 jars with respective lids and solvent resistant slide rack(s) with precise slide holding slots for 24 slides. The staining station is designed to hold the lids at the bottom of each jar to prevent mixing and/or losing them. The unique feature of this stand is that the user can add multiple staining stands to the existing one allowing the user to create their own stand as per the requirement. The linear staining station is compact and can fit in small areas. The durable material used in the frame gives an extra strength and sustains rough handling. It has been designed to save on reagent costs and helps in reducing the cost per slide.



Ordering Information	
<b>Description</b>	<b>Catalog#</b>
Deparaffinization & Staining Station (4 jars)	: SC004
Deparaffinization & Staining Station (8 jars)	: SC008

Technical Specifications for 4 Jars Staining Station	
Body Type	Plastic
Size (in cms) - Jars 4 capacity	33.5(L)×13.8(W)×14.1(H)
Weight - Jars 4 capacity	1.6kg~1.8kg / 3.5lb~3.9lb


4 Jars Staining Station Includes	
Staining Jars including Lids	: 4 Nos
Slide Rack	: 4 Nos
Stand (Upper Frame & Lower Frame)	: 2 Nos
Supporting Legs	: 4 Nos

# SLIDE STAINING CHAMBER



24 Slides Capacity

The multipurpose staining chamber is suitable for various slide staining protocols. The dark colour lid supports for light sensitive stains which are mainly used in immunohistochemistry, immunocytochemistry, immunofluorescence techniques etc. The chamber contains 4 raised rails to hold upto 24 slides and prevent the interference from the drain/waste from washings. The design allows to stack-up multiple racks to avoid space constrains.



Ordering Information	
<b>Description</b>	<b>Catalog#</b>
Slide Chamber	ST024-B

Technical Specifications	
Body Type	: Plastic
Size (in cms)	: 46×25×4
Slide Capacity	: 24

## CONTACT US

### USA

1257, Quarry Ln, Suite150  
Pleasanton, CA-94566, USA  
Phone: +1 925-218-6939

### Email

techsupport@pathnsitu.com  
customerservice@pathnsitu.com

### India

PS Heights, 2nd Floor, Plot# A-28/1/11B/B,  
Road#15 Nacharam IDA,  
Hyderabad-500 076, Telangana.  
Phone: 040-2701-5533 /5544  
Fax: 040-2701-5544