


# WATER BATH

## PS-WB



### Features

- Provides high and quick heat transmission
- Precise in temperature
- Programmable temperature settings
- Bright LED digital display
- Built-in battery backup system for memorized settings
- Black surface provides good visibility to identify the tissue sections
- Easy to clean and maintain

## Description

Water bath/Tissue Floatation Bath is an important step after cutting the paraffin sections and before the sections are placed on slides. It allows the tissue section to smoothen out before they are taken on the glass slide. Water bath ensures removal of any wrinkles and folds that are formed before placing them on the slide. The PS-WB Model water bath comprises of a surface that provides quick heat transmission. The temperature can be selected between the room temperature and 90 degrees.

## Technical Specifications

<b>Capacity</b>	1.5 litres
<b>Temperature Range</b>	Room temperature to 90°C ± 2°C
<b>Working Voltage</b>	AC220V / 50~60Hz
<b>Power</b>	200W
<b>Dimensions (in mm)</b>	350(L)x340(W)x120(H)
<b>Weight</b>	5kg

## Ordering Information

<b>Equipment Name</b>	<b>Catalog No</b>
Water Bath/Tissue Floatation Bath	PS-WB

## Contact Us

### USA

1257, Quarry Ln, Suite150  
Pleasanton, CA-94566, USA  
Phone: +1 925-218-6939

### Email

techsupport@pathnsitu.com  
customerservice@pathnsitu.com

[www.pathnsitu.com](http://www.pathnsitu.com)

### India

PS Heights, 2nd Floor, Plot# A-28/1/11B/B,  
Road#15 Nacharam IDA,  
Hyderabad-500 076, Telangana.

Phone: 040-2701-5533 /5544  
Fax: 040-2701-5544

